



VELNEZ®



CE
2460

NAME & ADDRESS OF THE MANUFACTURER:

Datt Mediprodukt Private Limited
Corporate Office: 56, Community Centre, East of Kailash, New Delhi-110065, India
Factory: Plot no: 52-54, 63 & 64 Roz-ka-Meo, Industrial Area Nuh, Distt. Mewat, Sohna, Haryana - 122103,
India www.dattmedi.com

EC REP

mdi Europa GmbH
Langenhagener
Str. 71
D-30855 Hannover-Langenhagen
Germany Phone: +49-511-3908 95312
Internet: www.mdi-europa.com
Trade Register Hanover HRB 58423

CATEGORY

Ear Pack

TRADE NAME

VelNez®

DESCRIPTION

VelNez® is a biodegradable composite that degrades within a few days on application and is intended as an ear dressing. VelNez® degrades within 20 days of application. It will also reduce fibrosis and at the same time promote healing and reduce the clotting time. During ear procedures, it can be applied to control compression of tissue, middle ear structure and prevent adhesion.

INTENDED USE

VelNez® is intended to be used as a packing dressing in the middle and external ear cavity after surgery. It acts as a space occupying dressing, preventing adhesion by separating the compromised mucosal surfaces. It also helps in wound healing and achieving haemostasis.

INSTRUCTIONS FOR USE OF VELNEZ® NASAL/ EAR DRESSING

Step 1: Removal from Package – Open the pack by removing the paperfilm and take out VelNez®.

Step 2: Pre-insertion – If required cut the VelNez® using sterile apparatus.

Step 3: Insertion – Using sterile apparatus, place the VelNez® inside ear cavity at the site of application.

Step 4: Confirm - Confirm that the device is placed at the intended position. VelNez® will fill itself by absorbing the surrounding fluids and blood.

CONSTITUENTS

Gelatin (BSE/TSE agent free) Chitosan, Polyvinyl Alcohol, Psyllium Husk

WARNING & PRECAUTIONS

- Do not re-sterilize.
- Do not re-use.
- Use the device prior to the "Use by (M)" specified on the package.
- To be used by Registered Medical Practitioner or a hospital
- Do not use if package is damaged.
- Discontinue the use if the application area shows signs of infection, irritation, maceration (whitening of surrounding skin) or itching & consult a healthcare professional.
- Should not be used for the patients with coagulation disorders.

RISK OF REUSE

VelNez® is supplied sterile for single use only. Do not re-sterilize as this may change product characteristics and lead to failure of intended use of the product. Damaged, unused, opened packages of VelNez® should be discarded as per disposal instructions. Sterility of the product is not guaranteed when the packaging is damaged or unused portion of the device is left in the pack

INDICATIONS

As ear pack for hemorrhage control or/and prevention of adhesion between the compromised surfaces after an ear surgery, traumatic bleeding, Tympanoplasty, Tympanostomy or Myringoplasty and other outer and middle ear surgeries.

HOW SUPPLIED

VelNez® is supplied as a sterile, single-use product in one pack.

CONTRAINDICATION

No known contraindications and adverse reactions reported.

Note: Do not use VelNez® if you are allergic (Hypersensitive) to its constituents.

NOTICE TO THE USER/ PATIENT

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority in which the user and/or patient is established

DISPOSAL OF DEVICE

Discard as per the hospital biomedical waste management system meeting local or national regulatory laws.

PATIENT POPULATION

Patients with age group from 18-60 years of age.

INTENDED USER

Surgeon, ENT specialist

STORAGE

VelNez® should be stored at a temperature between 4°C to 30°C. Keep away from direct sunlight. Store in a dry place.

SHELF LIFE

3 Years from the date of Mfg.

STERILIZATION METHOD

Gamma Irradiation

DEFINITION OF SYMBOLS USED

Symbol	Description
	Lot / Batch Number
	Manufacturing Date
	Use by
	Manufacturer
	Do not Reuse
	Do not Re-sterilize
	Sterilized using Irradiation (Gamma Radiation)
	Do not use if package is damaged
	Keep away from sunlight
	Store in dry place
	Storage Temperature
	Caution
	Instructions for use
	Latex Free
	Catalogue Number
	Animal origin
	Medical Device